

# Isotherm.

# Low Temperature Incubators

Introducing Esco Isotherm<sub>®</sub> - world class Low Temperature Incubators from Esco with the technologies and compliance to prove it. Ergonomic, intuitive interfaces, microprocessor PID controls with programming options, pre-heat chamber technology, dual auto-defrosting system, UV disinfection, precisely tuned and tested ventilation and insulation package, all supported by Esco's solutions - based sales and service representatives worldwide.







### Pre-heat chamber technology

No exposed heating elements located inside the chamber to ensure maximum user safety.

> Stable heating and maximum temperature uniformity in the chamber.

Standard temp setting range 0°C up to 100°C for maximum application flexibility.

Secure 2-point door seal and eccentric hinge ensure maximum gasket compression for stable chamber temperature.







#### Ventilation system

Forced convection design allows rapid temperature response rates, improves uniformity, and reduces fluctuation.

Fan speed and air exchange rates are adjustable. Low noise during operation.

Ventilated stainless steel shelves contribute to uniform air circulation..


#### **Auto-defrosting system**

Auto-heating activates and continues for a predetermined time during opration.

Auto-defrosting during operation

Auto-defrosting activates regularly. Influence on temperature fluctuation and uniformity is minimal.





#### **Easy To Service**

Diagnostic functions in the microprocessor include historical read-out of temperatures.

Diagnostic menu provides read-out of all sensor inputs and controller settings.

Service can be carried out from the front.

All electronics components are isolated from the work chamber and easily accessible for replacement.

Low service costs.





#### **Side Access port**

For temperature validation & mapping.

Isotherm.



#### **Microprocessor PID Control Technology**

Soft touch controls for temprature and UV are easy to clean. Tuned PID control ensures fast ramp time, prevents over-shoot and ensures stable temperature once setpoint is achieved.




#### **UV** disinfection

Auto-running disinfection after startup and during testing process as well, is adjustable to meet different requirements. UV light function stops when door is opened and resumes after it is closed.

# **Key Features**

- Pre-heat chamber technology.
  - No exposed heating elements located inside the chamber to ensure maximum user safety.
  - Stable heating and maximum temperature uniformity in the chamber.
  - Standard temp setting range 0° C up to 100° C for maximum application flexibility.
  - Secure 2-point door seal and eccentric hinge ensure maximum gasket compression for stable chamber temperature.
- Ventilation system.
  - Forced convection design allows rapid temperature response rates, improves uniformity, and reduces fluctuation.
  - Low noise during operation.
  - Ventilated stainless steel shelves contribute to uniform air circulation.
- Microprocessor PID Control Technology.
- Auto-defrosting system
  - Auto-heating activates and continues for a predetermined time during opration.
  - Auto-defrosting during operation -Auto-defrosting activates regularly.
    Influence on temperature fluctuation and uniformity is minimal.
- UV disinfection
  - Auto-running disinfection after startup and during testing process as well, is adjustable to meet different requirements.
  - UV light function stops when door is opened and resumes after it is closed.
- Easy-to-Clean.
- Easy-to-Service.



#### **Built-In Protection**

External surfaces are powder coated with Esco **ISOCIDE**" to eliminate 99.9% of surface bacteria within 24 hours of exposure.



# General Specifications, Isotherm Low Temperature Incubators (IFC)

Model			IFC-110-8	IFC-240-8
Volume			110 litre (3.88 cu.ft)	240 litre (8.48 cu.ft)
Temperature Range			0°C ~ 100°C	
Temperature Variation	At 10°C(±K)		≤±0.6°C	≤±0.6°C
	At 15°C(±K)		≤±0.4°C	≤±0.4°C
	At 25°C(±K)		≤±0.3°C	≤±0.3°C
	At 37°C(±K)		≤±0.3°C	≤±0.5°C
Temperature Fluctuation	At 10°C(±K)		≤±0.3°C	≤±0.3°C
	At 15°C(±K)		≤±0.3°C	≤±0.3°C
	At 25°C(±K)		≤±0.3°C	≤±0.3°C
	At 37°C(±K)		≤±0.3°C	≤±0.3°C
Heating up time to 37°C from the ambient temperature			31 minutes	37 minutes
Heating up time to 50°C from the ambient temperature			33 minutes	50 minutes
Recovery time after door was opened for 30 sec	At 5°C		3 minutes	5 minutes
	At 37°C		2 minutes	3 minutes
	At 50°C		2 minutes	3 minutes
Cooling down time from 22°C to 0°C			90 minutes	120 minutes
Cooling down time from 22°C to 5°C			60 minutes	80 minutes
Cooling down time from 22°C to 10°C			34 minutes	48 minutes
Heat emission at 37°C set point (compressor on) (Watt)			217	238
Heat emission at 37°C set point (compressor off) (Watt)			61.8	80.7
Electrical (200-240V, AC, 50/60Hz, 1Ø)	Power*	at 15°C	400 W	481 W
		at 25°C	431 W	563 W
	Cabinet Full Load Amps (FLA)		6 A	6 A
	Cabinet BTU	Set Point 15°C	1364.84	1641.23
		Set Point 25°C	1470.63	1921.03
Incubator Construction	Main Body		Electro galvanized steel with white oven baked epoxy powder-coated finish	
	Chamber		Stainless steel, grade 304	
Number of Shelves	Standard		2	2
	Maximum		4	8
Load Per Shelf		30 kg (13.6 lbs)		
Max. Total Load			60 kg (27.3 lbs)	
External Dimensions (W x D x H)			820 x 730 x 1185 mm (32.3" x 28.7" x 45.6")	841 x 871 x 1462 mm (33.1" x 34.3" x 53.3")
Internal Dimensions (W x D	x H)		600 × 480 × 399 mm (23.6" x 18.9" x 15.7")	645 x 700 x 530 mm (25.4" x 27.6" x 20.9")
Net Weight			134 Kg (295.41 lbs)	164 kg (361.55 lbs)
Shipping Weight			166 Kg (365.96 lbs)	195 kg (429.90 lbs)
Shipping Dimensions, Maximum (W x D x H)			878 x 787 x 1425 mm (34.5" x 30.9" x 56.1")	891 x 933 x 1628 mm (35" x 36.7" x 64.1")
Shipping Volume, Maximum			0.98 m³ (34.60 cu.ft)	1.35 m <sup>3</sup> (47.67 cu.ft)

- All technical specifications are specified for units with standard equipment at an ambient temperature of 25°C and a voltage fluctuation of  $\pm 10\%$ . The temperature data is determined in accordance to DIN 12880 standards. All indications are average values, typical for units produced in series.
- Esco reserves the right to alter technical specifications at all times.
- In order to calculate the current at maximum power consumption, divide maximum power consumption by the voltage.





Esco Technologies, Inc. • 2940 Turnpike Drive, Units 15-16 • Hatboro, PA 19040, USA Toll-Free USA and Canada 877-479-3726 • Tel 215-441-9661 • Fax 215-441-9660 us.escoglobal.com • usa@escoglobal.com

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777 Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escoglobal.com www.escoglobal.com

Esco Global Offices | Beijing, China | Kuala Lumpur, Malaysia | Manama, Bahrain | Guangzhou, China | Hanoi, Vietnam | Melaka, Malaysia | Mumbai, India | Philadelphia, PA, USA | Salisbury, UK | Shanghai, China | Seoul, Korea | Delhi, India | Osaka, Japan | Manila, Philippines | Midrand, South Africa | Jakarta, Indonesia | Singapore

